

## Advanced Lift Therapy CU-X

Ref.	4327 / 4327P
Description	Vitamin A cure with Copper Peptides
Retail size	Syringe with refills 2 x 10 ml, cure of app. 8 weeks
Professional size	30 ml airless dispenser
Skin type	For mature skin with wrinkles and sagging contours
Product galenics	Apricot colored serum



### DESCRIPTION

The latest breakthrough in anti-aging research.

Copper plays a key role in generating elastin, collagen and the skin's own hyaluronic acid, as well as in interweaving collagen and elastin fibres. On top of this, it also has anti-inflammatory properties. A lack of copper leads to a breakdown of skin regeneration. Inflammation can spread across the skin. This can lead to premature aging of the skin (inflamm-aging). This is where the inspira:med Copper Peptides help - special peptides work like a kind of cosmetic drone and bring the copper directly to the skin cells that urgently need it, in order to fight the inflammatory process effectively and regenerate the skin starting from the inside and working outwards.

inspira:med combines the effects of Copper Peptides with vitamin A, which stimulates skin regeneration. A highly efficient lifting peptide helps to visibly lift facial contours. In addition the skin is also supplied with moisture thanks to long and short-chain hyaluronic acids. Fine lines and wrinkles are softened.

The result: visibly younger and smoother skin with redefined contours.

### Effects and characteristics

- Combination of Copper Peptides and Vitamin A
- Stimulates generating elastin, collagen and the skin's own hyaluronic acid, as well as in interweaving collagen and elastin fibres
- Special peptides carry the copper into skin cells where it is needed
- Effectively fights the inflammatory process
- Regenerates the skin starting from the inside and working outwards
- Highly efficient lifting peptide
- Moisturizing
- Visible lifting of facial contours
- Visibly younger and smoother skin

### Special ingredients:

- **Copper Peptides:** Special peptides, work like a kind of cosmetic drone and bring the copper directly to the skin cells that urgently need it, in order to fight the inflammatory process effectively and regenerate the skin starting from the inside and working outwards.
- **Highly efficient lifting peptide:** Visibly lifts facial contours.
- **Long- and short-chain hyaluronic acid:** Are biotechnologically produced and have a fill-in effect and moisturize. Fine lines and wrinkles are softened.
- **Vitamin A:** Retinol, smoothes the skin. Increases the cell regeneration and skin elasticity. In the skin, vitamin A is converted into vitamin A acid, dissolves the top skin layer and stimulates the production of new skin cells.



# inspira:med

## HOME CARE APPLICATION

Before the first application: Unscrew the syringe at the top, remove the cap of the refill cartridge, place the cartridge in the syringe and screw back up, unlock the syringe, pull the cap off the top and draw up the product. To be used as a cure for up to 8 weeks, apply in the evening after cleansing and before applying overnight products.

## PROFESSIONAL APPLICATION

Apply 2-3 drops onto the face, neck and neckline and allow working in briefly.

Then apply one of the following masks: Calming Gel Mask, Bio Cellulose Mask Age Defying or Peel Off Mask Filmalgine®.

